



# PATHWAY<sup>®</sup> STM-10

VAGINAL AND RECTAL INTRACAVITY STIMULATOR

Regaining Pelvic Floor Function through Sensory Awareness



# PATHWAY<sup>®</sup> STM-10



## FEATURES:

- Transition from Clinical to Home Treatment
- One Channel of Vaginal or Rectal Intracavity Stimulation
- 12.5, 50, Sequential 12.5 / 50, 100 and 200Hz Frequencies
- Microprocessor Controlled
- Powered by a Readily Available 9V Battery
- Programmable Treatment Time and On/Off Intervals
- Treatment Protocol Storage
- Treatment Compliance Monitor
- Emulates the Stimulation Capabilities of its Clinical Counterparts

## TECHNICAL SPECIFICATIONS:

**Pulse Amplitude (mA):** 1-99

**Frequency Rate (Hz):** 12.5 / 50 / 100 / 200

**Duty Cycle (on/off) (S):** 5/5, 5/10, 10/10, 10/20, Continuous

**Session Time:** 5, 10, 15, 20, 25, 30 Minutes

**Ramp Up (S):** 2 Fixed

**Ramp Down (S):** 1 Fixed

**Pulse Width (uS):** 300

**Charge Per Pulse (uC):** 17

**Power Source:** 9V Alkaline Battery

**Dimensions (In):** 4.75 x 2.5 x 1

**Weight (oz):** 4.2 (5.8 with 9V Alkaline Battery)

**Compliance Use Timer:** Hours:Minutes (99:59 Max)



ISO 13485 CERTIFIED • © The Prometheus Group®, April 2016. All rights reserved.



Scan here to learn what The Prometheus Group® can do for you.

Contact us today for more information or to schedule a demonstration  
info@theprogrp.com | www.theprogrp.com | 1.800.442.2325